

App Guide  
Ground Reaction Force(GRF)  
Bilateral

SageMotion  
Wearable Biofeedback System





## **Table of Contents**

Components

Wirelessly Connect to Computer or Cellphone

GRF Bilateral APP

# Components



Hub



Nodes (8x)



Battery



Node Straps: *Medium (8x), Short (4x), Long (2x)*



Cable A (10x)

-Connect Hub to Battery  
-Charge Nodes & Battery



Cable B (*optional use*)

-Connect Hub to Computer




Node Charging Station

# Wirelessly Connect to Computer or Cellphone

## 1) Connect Cable A to Battery and to Hub




## 2) On Computer/Cellphone, Connect to Wi-Fi: "SageMotion"



*Note 1: Need to wait for up to 1 minute for "SageMotion" to appear in Wi-Fi list. If it doesn't appear, try turning the Wi-Fi off and then on again on the computer/cellphone.*

*Note 2: Hub is connected after clicking "Connect" even if in Windows it shows "Connecting" or "No internet, open".*

## 3) On Computer/Cellphone, in Chrome Address Bar, Go To <http://192.168.12.1>



**[Note] If Computer Doesn't Have Wi-Fi:** plug in Cable B to the Hub and to the ethernet port of your computer, then in chrome address bar, go to **<http://192.168.137.1>**

# GRF Bilateral App

The purpose of this GRF Bilateral app is to predict and stream Ground Reaction Forces (GRFs) in real-time based on IMUs data, facilitating the analysis of force distribution during various physical activities.

## 1) Turn on 7 Nodes




*Slide switch toward middle to turn node on*



*Green light will blink after the node is on and running*

## 2) Select “GRF Bilateral” App



## 3) Click “Search”

Node List



# GRF Bilateral App (cont.)

4) Configure 7 Sensor Nodes as Shown Below:

Node List ?



Type	Position	MAC					
sensor	L FOOT	18:04:ED:E7:83:4F			Blink	Vibrate	+
sensor	R FOOT	18:04:ED:EB:48:96			Blink	Vibrate	+
sensor	R SHANK	18:04:ED:EB:21:A1			Blink	Vibrate	+
sensor	R THIGH	18:04:ED:F9:A1:7D			Blink	Vibrate	+
sensor	PELVIS	18:04:ED:EB:2B:D8			Blink	Vibrate	+
sensor	L SHANK	18:04:ED:7E:65:C6			Blink	Vibrate	+
sensor	L THIGH	18:04:ED:EB:23:33			Blink	Vibrate	+

Number of Nodes Required = 7

5) Click “Connect Nodes”



6) “Ready to collect data” Will Appear after Node Connection is Complete



GRF\_Online



Start

Stop




Ready to collect data



# GRF Bilateral App (cont.)



7) Attach the foot nodes by sticking the rough side of the Velcro to the shoes and the soft side of the Velcro to the nodes. Then, fasten the Velcro together



*For both nodes, the on/off switch points to the body*

8) Thread Straps through the Remaining 5 Nodes

## How to Thread Straps






# GRF Bilateral App (cont.)

9) Attach Nodes to Anterior Shanks – Midway between Femur Lateral Epicondyle and Fibula Apex of Lateral Malleolus

*For both nodes, the on/off switch points upwards*




*Femur Lateral Epicondyle*

*Fibula Apex of Lateral Malleolus*

10) Attach Nodes to Anterior Thighs – Midway between Femur Greater Trochanter and Femur Lateral Epicondyle

*For both nodes, the on/off switch points upwards*



*Femur Greater Trochanter*

*Femur Lateral Epicondyle*

# GRF Bilateral App (cont.)

11) Attach Nodes to Anterior Pelvis – Midway between Left and Right Anterior Superior Iliac Spine

*The on/off switch points upwards*



*Anterior Superior Iliac Spine*

12) Click “Blink” for each Node to Confirm Correct Locations (red LED for given node blinks 3 times on click)

## Node List ?



Type	Position	MAC					
sensor	L FOOT	18:04:ED:E7:83:4F			Blink	Vibrate	
sensor	R FOOT	18:04:ED:EB:48:96			Blink	Vibrate	
sensor	R SHANK	18:04:ED:EB:21:A1			Blink	Vibrate	
sensor	R THIGH	18:04:ED:F9:A1:7D			Blink	Vibrate	
sensor	PELVIS	18:04:ED:EB:2B:D8			Blink	Vibrate	
sensor	L SHANK	18:04:ED:7E:65:C6			Blink	Vibrate	
sensor	L THIGH	18:04:ED:EB:23:33			Blink	Vibrate	

Number of Nodes Required = 7

## GRF Bilateral App (*cont.*)

13) Enter App Configuration Settings including the Subject's Height and Weight

### App Configuration

Trial Name	<input type="text" value="trial_3"/>
<b>Save Options</b>	
Save Mode	<input type="text" value="csv"/> ▼
<b>User Info</b>	
Subject Weight (Kg)	<input type="text" value="59"/>
Subject Height (m)	<input type="text" value="1.75"/>

# GRF Bilateral App (cont.)

14) Click “Start” to Start Running the App



15) After the Trial is Finished, Click “Stop”



16) After Clicking “Stop”, a File from that Trial will appear under Data Management. Click the file (e.g.trial\_3) to download it to the device

## Data List

<input type="checkbox"/> Name	Date▲	Duration	App			
<input type="checkbox"/> <u>GRF Bilateral trial 3</u>	2024-07-24-15-34-40	0:00:16	GRF Bilateral			

# GRF Bilateral App (*cont.*)

## Description of Data in Downloaded File

**time (sec):** time since trial start

**plate\_1\_force\_x:** X component of the GRF acting on the left foot (horizontal force in the medial-lateral direction).

**plate\_1\_force\_y:** Y component of the GRF acting on the left foot (horizontal force in the anterior-posterior direction).

**plate\_1\_force\_z:** Z component of the GRF acting on the left foot (vertical force in the superior-inferior direction).

**plate\_2\_force\_x:** X component of the GRF acting on the right foot (horizontal force in the medial-lateral direction).

**plate\_2\_force\_y:** Y component of the GRF acting on the right foot (horizontal force in the anterior-posterior direction).

**plate\_2\_force\_z:** Z component of the GRF acting on the right foot (vertical force in the superior-inferior direction).

**SensorIndex\_1/2/3/4/5/6/7:** index of raw sensor data

**AccelX/Y/Z\_1/2/3/4/5/6/7** (m/s<sup>2</sup>): raw acceleration data

**GyroX/Y/Z\_1/2/3/4/5/6/7** (deg/s): raw gyroscope data

**MagX/Y/Z\_1/2/3/4/5/6/7** (μT): raw magnetometer data

**Quat1/2/3/4\_1/2/3/4/5/6/7:** quaternion data

**Sampletime\_1/2/3/4/5/6/7:** timestamp of each sensor

**Package\_1/2/3/4/5/6/7:** package number of each sensor